



# Gemlogis VISTA



The Gemlogis VISTA is a portable diamond verification instrument for loose and mounted diamonds, designed to separate diamonds from synthetic diamonds. The VISTA is manually operated and gives an automatic test result.

The VISTA automatically makes a measurement and displays one of the following results:

- ◇ Natural Diamond
- ◇ Type II Diamond

## Stone Testing Capabilities

Weight of stones:	<b>0.02 ct - 12 ct</b>
Size:	<b>1.5 mm - 15.3 mm</b>
Colour of stones:	<b>D-J</b>
Shape of stones:	<b>All</b>
Diamond simulants:	<b>No</b>
Mounted stones:	<b>Yes</b> closed and open back mountings, not yet ASSURE Tested

## Instrument Capabilities

Automated feed:	<b>No</b>
Automated results:	<b>Yes</b>
Automated dispense:	<b>No</b>
Detect or Refer (synthetic diamonds)	<b>Refer</b>
Detect or Refer (diamond simulants)	<b>N/A</b>
Multiple stones at once:	<b>No</b>
Training required:	<b>No</b>

## Instrument Performance

The results demonstrate the performance of the Gemlogis VISTA tested by a novice operator:

	<b>Core Sample</b>
<b>Diamond False Positive Rate</b> <i>Optimal rate 0%</i>	<b>0%</b>
<b>Diamond Referral Rate</b> <i>Optimal rate 0%</i>	<b>7.5%</b>
<b>Diamond Accuracy Rate</b> <i>Optimal rate 100%</i>	<b>92.5%</b>
<b>Speed</b> <i>Stones tested per hour</i>	<b>212</b>

## Instrument Specifications

 **4.1 cm (w) x 9.8 cm (d) x 11 cm (h)**

 **0.3 kg**

**599 USD** August 2019

[gemlogisusa.com/gemlogis-vista.html](http://gemlogisusa.com/gemlogis-vista.html)

[dana@gemlogisusa.com](mailto:dana@gemlogisusa.com)

**+1 800 833 1090**

**ASSURE Tested September 2019**

**ASSURE ID 94513**